

Ring Pull Strain Analysis (RPSA) Tutorial

Introduction


This tutorial is intended for a novice user to be able to run the Ring Pull Strain Analysis program.

Getting Started

Please download Python from anaconda.org. There are multiple ways to download Python, however, this one contains all the necessary modules. Features of Anaconda will be displayed heavily in this tutorial.

File Structure

The files that are needed to run this program are shown below




Below are some images showing the type of data present in each file. Note that the .csv files must have specific headers.

Load Frame File:

Each data point has to have image files associated with it. Specify the exact names of the image here. The default/first images should go in the top_img_file column. The side_img_file column is for coating analysis. The filenames must have the extension listed as well. You may either specify the FULL filepath OR you may have all the images in the same folder as this file.

	A	B	C	D	E
1	top_img_file	side_img_file	time (sec)	displacement (mm)	load (N)
2	img0-08312022144721-9.tiff	img1-08312022144721-7.tiff	0	0	1.52449
3	img0-08312022144722-10.tiff	img1-08312022144722-8.tiff	1	0.005853	3.217074
4	img0-08312022144723-11.tiff	img1-08312022144723-9.tiff	2	0.015224	3.510204
5	img0-08312022144724-12.tiff	img1-08312022144724-10.tiff	3	0.02253	3.897891
6	img0-08312022144725-13.tiff	img1-08312022144725-11.tiff	4	0.030023	4.361225
7	img0-08312022144726-14.tiff	img1-08312022144726-12.tiff	5	0.037377	4.739456
8	img0-08312022144727-15.tiff	img1-08312022144727-13.tiff	6	0.04475	5.20279

Image Files:



csv files:


- These files are most often generated from a third party digital image correlation software. Supported file modes are VIC-2D 6, VIC-2D 7, and DICEngine.

Necessary columns


A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R	S	T	U	V	W	X	Y
"x"	"y"	"u_c"	"v_c"	"u"	"v"	"u_c"	"v_c"	"exx"	"eyy"	"exy"	"fxx"	"fyy"	"fxy"	"fxx"	"fyy"	"e1"	"e2"	"gamma"	"e_tres"	"e_vonmi"	"sigma"	"Xp"	"Yp"	"Zp"
2	24	322	-1173	702	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	-1	0	0	0
3	55	322	-1169	702	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	-1	0	0	0
4	59	322	-1165	702	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	-1	0	0	0
5	63	322	-1161	702	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	-1	0	0	0
6	67	322	-1157	702	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	-1	0	0	0
7	71	322	-1153	702	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	-1	0	0	0
8	75	322	-1149	702	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	-1	0	0	0
9	79	322	-1145	702	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	-1	0	0	0

Opening and running the script


Start by opening Spyder:



It is required to change the graphics backend to display the plots in a separate window. The backend cannot be inline, so any other setting should work. This will allow Spyder to create additional windows for any plots or figures that are made.



Now the RPSA example script can be opened.



Now, walking through the script. Here are some imported packages. If you downloaded Anaconda, these should all be available with minimal hassle. We also define pi as we will often use it. Note the red box – this is where we import the module RingPullStrainAnalysis.py from a different folder.

```
9  ## import the module from a different folder
10 import sys
11 sys.path.insert(0, 'E:\\Projects\\RingPuLL\\RPSA')
12 from RingPullStrainAnalysis import *
13
14 import numpy as np
15 pi = np.pi
16 import matplotlib.pyplot as plt
17
```

Here are the variables you can change in the code. The most important are the directory and the filename for the load frame data. Here is also where you specify geometric dimensions of the ring.

This variable must have the full filepath	Requires either double backslashes (shown) or single frontslash
<pre>21 22 ## set all the variables you will need to run this RingPullStrainAnalysis code: 23 24 ## the file with the load frame data and their corresponding images 25 LF_file = 'E:\\Projects\\RingPull\\RPSA_example\\sample_data\\sample_data_LF.csv' 26 #The DIC analysis software that was used 27 DIC_software = 'VIC-2D' 28 ## geometric dimensions of the ring pull test 29 d_mandrel = 3.0 30 OD = 10.3 31 ID = 9.46 32 W = 1</pre>	

All filenames must have their extensions on them


This line of code creates an instance of the RingPull class which is intended to encompass an entire test. It loads in the load frame data and stores it in memory.

```
37
38  ## create RingPull object
39  test = RingPull(LF_file=LF_file,
40                  DIC_software=DIC_software,
41                  ID=ID, OD=OD, d_mandrel=d_mandrel, W=W,
42                  get_geometry_flag=True)
43
```

You can put a digital extensometer on the sample to get an adjusted displacement. This pops up a user interface allowing the user to select two points to track. This will put an additional column in the load-displacement DataFrame internally in Python. The user can save this data externally to a .csv file for future use.


```
45
46  ## use the DIC data to create a digital extensometer on the test sample
47  test.digital_extensometer()
48
49  ## you can also save the adjusted displacement calculation for future use
50  test.save_data()
51
52  ## The complete load-displacement data with calculated stress and strain values
53  ## can be called by looking at the state variable, df
54  df = test.df
55
```

Choose 2 points to create a digital extensometer on the sample
Right click to finish.
Middle mouse to remove points




Strain values at arbitrary points can be probed. Here is an example of strain values being obtained from the DIC mesh and plotted.

```
57  
58  ## plot the strain distribution from one of the DIC images  
59  n = 245  
60  theta = pi/2+1e-4  
61  a = np.linspace(0,1,50)  
62  e = test.open_Image(n).get_value(a,theta,mode='ett', extrap=False)  
63  
64  f,ax = make_figure()  
65  ax.plot(a,e)  
66  ax.set_xlabel('location on ring (a)')  
67  ax.set_ylabel('strain')  
68
```




You can also analyze this curve as if it were a tensile curve. When this command is run, user graphical input is requested.

```
68
69  ## analyze the curve as if it were a stress strain curve from a tensile test
70  ## and output important material parameters
71  E,YS,UTS,eps_u,eps_nu,eps_total,toughness = test.process_stress_strain_curve()
72
```




```
73
74  ## Some more plotting methods
75  ax = test.plot_stress_strain()
76  ax.set_xlabel('Strain [mm/mm]')
77  ax.set_ylabel('Stress [MPa]')
78
```



Now, we can open a single image file to look at a specific point in the test. You can also access the data in spyder's variable editor.

We opened a random image (number 354). Note that the image numbers correspond with their place in the load frame file, and NOT with any numbers in their filenames

```
60
61  ## open one of the DIC_image classes from the RingPull object
62  img = test.open_DIC(354)
63  print(type(img)==DIC_image)
64
65  ## again, this data is pulled in from the csv file. You can see the data here:
66  df2 = img.df
67
```

Now we can also plot a few more things, shown below.

```
89
90  ## Plots the DIC results overlayed on the image
91  img.plot_Image(state='reference')
92  img.plot_Image(state='deformed')
93
94  ## You can also plot strains in polar coordinates.
95  img.plot_Image(state='deformed',mode='ett')
96
97
```

